

# The Leeds Inhaler Device Guide: Inhaler Technique Instructions for Healthcare Professionals and Patients

1st Edition  
December 2015

**Authors:** Dr Toby Capstick, Lead Respiratory Pharmacist  
Kelly Atack, Specialist Clinical Pharmacist

**Reviewers:**

## Contents

How to Use the Leeds Inhaler Device Guide	3
Types of Inhaler Device	3
Use of placebo inhalers	4
Use of In-Check DIAL Inspiratory Flow Meters	4
Other Tools Available to Assess Inspiratory Flow	4
Recommendations for Performing Inhaler Technique Assessments	4
Aerosol science & inhaled medication	5
Inhaler Device Monographs:	
• Accuhaler	7
• Aeroliser	8
• Autohaler	9
• Breezhaler	10
• Clickhaler	11
• Easi-Breathe	12
• Easyhaler	13
• Ellipta	14
• Forspiro	15
• Genuair	16
• HandiHaler	17
• NEXThaler	18
• Novolizer	19
• pMDI (Pressurised Metered Dose Inhaler)	20
• pMDI + Spacer (Multiple Breath Method)	21
• pMDI + Spacer (Single Breath Method)	22
• Podhaler	23
• Pulvinal	24
• Respimat	25
• Spiromax	26
• Turbohaler	27
• Turbospin	28
• Twisthaler	29
Monograph References	30

## Glossary of Drug Classes

SABA	Short-Acting Beta <sub>2</sub> -Agonist
SAMA	Short-Acting Muscarinic Antagonist
LABA	Long-Acting Beta <sub>2</sub> -Agonist
LAMA	Long-Acting Muscarinic Antagonist
ICS	Inhaled Corticosteroid

The following brand names and registered trademarks mentioned throughout this guide are acknowledged:

### Inhaler Devices:

Accuhaler<sup>®</sup>, Aeroliser<sup>®</sup>, Autohaler<sup>®</sup>, Breezhaler<sup>®</sup>, Clickhaler<sup>®</sup>, Easi-Breathe<sup>®</sup>, Easyhaler<sup>®</sup>, Ellipta<sup>®</sup>, Forspiro<sup>®</sup>, Genuair<sup>®</sup>, HandiHaler<sup>®</sup>, NEXThaler<sup>®</sup>, Novolizer<sup>®</sup>, pMDI, Podhaler<sup>®</sup>, Pulvinal<sup>®</sup>, Respimat<sup>®</sup>, Spiromax<sup>®</sup>, Turbohaler<sup>®</sup>, Turbospin<sup>®</sup>, Twisthaler<sup>®</sup>, AeroChamber Plus<sup>®</sup>, AeroChamber Plus for Infant or Toddler<sup>®</sup>, Able Spacer<sup>®</sup>, A2A Spacer<sup>®</sup>, DispozABLE Spacer<sup>®</sup>, Optichamber Advantage<sup>®</sup>, Optichamber Diamond<sup>®</sup>, Pocket Chamber<sup>®</sup>, Space Chamber Plus<sup>®</sup>, Volumatic Spacer<sup>®</sup>, Vortex<sup>®</sup>

### Brand names:

AirFluSal<sup>®</sup>, Airomir<sup>®</sup>, AirSalb<sup>®</sup>, Anoro<sup>®</sup>, Asmabec<sup>®</sup>, Asmanex<sup>®</sup>, Asmasal<sup>®</sup>, Atimos Modulite<sup>®</sup>, Atrovent<sup>®</sup>, Bricanyl<sup>®</sup>, Budelin<sup>®</sup>, Clenil Modulite<sup>®</sup>, Colobreathe<sup>®</sup>, Duaklir<sup>®</sup>, DuoResp<sup>®</sup>, Eklira<sup>®</sup>, Flixotide<sup>®</sup>, Flutiform<sup>®</sup>, Foradil<sup>®</sup>, Fostair<sup>®</sup>, Incruse<sup>®</sup>, Neoven<sup>®</sup>, Onbrez<sup>®</sup>, Oxis<sup>®</sup>, Pulmicort<sup>®</sup>, Qvar<sup>®</sup>, Relvar<sup>®</sup>, Salamol<sup>®</sup>, Salbulin<sup>®</sup>, Seebri<sup>®</sup>, Seretide<sup>®</sup>, Serevent<sup>®</sup>, Sirdupla™, Spiolto<sup>®</sup>, Spiriva<sup>®</sup>, Striverdi<sup>®</sup>, Symbicort<sup>®</sup>, Tobi<sup>®</sup>, Ultibro<sup>®</sup>, Ventolin<sup>®</sup>, Vertine<sup>®</sup>.

### Photo credits

Acknowledgement for the use of promotional product images is specified in each monograph. Where no promotional product images were available, photos of placebo inhalers are used.

## How to Use the Leeds Inhaler Device Guide

Prior to prescribing any new inhaler, each patient should receive training and education in the use of the inhaler device and have demonstrated satisfactory inhaler technique.<sup>1,2,3</sup> Additionally, correct inhaler technique should be assessed regularly at each clinic visit, and following an exacerbation.<sup>1,2</sup> If a patient is unable to use a particular device correctly, an alternative device should be used.<sup>1</sup>

Decisions about which inhaler device should be prescribed for any patient is complicated by the wide range of inhaler devices available in the UK market. As of the end of 2015, there are 21 different inhaler devices licensed for treating patients with asthma, COPD or cystic fibrosis. This guide aims to provide both healthcare professionals and patients with information about how to use each inhaler device, as well as describing key features of each device and essential 'need to know' details to explain the critical steps required to use each device correctly.

The Leeds Inhaler Device Guide does not aim to make any recommendation about which inhaler device should be prescribed for any patient. The choice of device to prescribe will depend on a number of factors, including indication, local formularies, how easy the inhaler device is for any individual to use, patient preference, and cost. However, it is important to remember that **the best inhaler device is the one that the individual patient can use**.

## Components of Inhaler Device Monographs

<table border="1"> <tr> <td>Type of inhaler:</td> <td></td> </tr> <tr> <td>In use expiry:</td> <td></td> </tr> <tr> <td>Available as</td> <td></td> </tr> <tr> <td>SABA</td> <td></td> </tr> <tr> <td>SAMA</td> <td></td> </tr> <tr> <td>LABA</td> <td></td> </tr> <tr> <td>LAMA</td> <td></td> </tr> <tr> <td>LAMA/LABA</td> <td></td> </tr> <tr> <td>ICS</td> <td></td> </tr> <tr> <td>ICS/LABA</td> <td></td> </tr> </table>	Type of inhaler:		In use expiry:		Available as		SABA		SAMA		LABA		LAMA		LAMA/LABA		ICS		ICS/LABA		<p>This section provides background information on the type of inhaler device in each monograph; whether it is a pressurised metered dose inhaler (pMDI), a breath-actuated pMDI, Soft Mist Inhaler, or a dry powder inhaler.</p> <p>If patients require more than one inhaler device to manage their respiratory disease, a list of drugs available in each device is provided, divided by class of drug. This may aid prescribing decisions, to ensure consistency of prescribing patterns.</p> <p>As many inhaler devices may have relatively short in-use life span, this information is provided in the 'Expiry' section.</p>
Type of inhaler:																					
In use expiry:																					
Available as																					
SABA																					
SAMA																					
LABA																					
LAMA																					
LAMA/LABA																					
ICS																					
ICS/LABA																					
<table border="1"> <tr> <td colspan="2">Instructions for Use</td> </tr> <tr> <td>1. Preparation</td> <td></td> </tr> <tr> <td>2. Priming</td> <td></td> </tr> <tr> <td>3. Exhaling</td> <td></td> </tr> <tr> <td>4. Mouth</td> <td></td> </tr> <tr> <td>5. Inhalation</td> <td></td> </tr> <tr> <td>6. Breath holding</td> <td></td> </tr> <tr> <td>7. Closing &amp; Repeating</td> <td></td> </tr> </table>	Instructions for Use		1. Preparation		2. Priming		3. Exhaling		4. Mouth		5. Inhalation		6. Breath holding		7. Closing & Repeating		<p>The Instructions for Use section of the monograph provides detailed guidance on exactly how to use each inhaler device.</p> <p>To provide a consistent approach to using inhaler devices, this has been standardised into 7 steps.</p> <p>The important critical steps that are required to obtain the most benefit from each inhaler device are highlighted in bold text.</p>				
Instructions for Use																					
1. Preparation																					
2. Priming																					
3. Exhaling																					
4. Mouth																					
5. Inhalation																					
6. Breath holding																					
7. Closing & Repeating																					
<table border="1"> <tr> <td colspan="2">Key Features</td> </tr> <tr> <td>Dexterity</td> <td>visually required</td> </tr> <tr> <td>Feedback: Dose Counter</td> <td></td> </tr> <tr> <td>Feedback: Tactile/Puff</td> <td></td> </tr> <tr> <td>Feedback: Sound</td> <td></td> </tr> <tr> <td>Feedback: Visual</td> <td></td> </tr> <tr> <td>Dose Lock-Out</td> <td></td> </tr> <tr> <td>Dose Consistency</td> <td></td> </tr> <tr> <td>Device Resistance</td> <td>officer resistance</td> </tr> <tr> <td>Inspiratory Flow</td> <td></td> </tr> </table>	Key Features		Dexterity	visually required	Feedback: Dose Counter		Feedback: Tactile/Puff		Feedback: Sound		Feedback: Visual		Dose Lock-Out		Dose Consistency		Device Resistance	officer resistance	Inspiratory Flow		<p>The Key Features section outlines different aspects of inhaler devices that may help healthcare professionals and patients understand how easy or hard different devices are to use.</p> <p>A subjective assessment has been made of the complexity of the inhaler device and the level of <i>Dexterity</i> that patients may need to have to operate it effectively. Devices described as requiring a high level of dexterity may not be suitable for people with physical impairment.</p> <p>Description of the <i>Feedback</i> mechanisms of each inhaler device illustrates ways in which patients may know that they have used their inhaler correctly, or know whether it is empty. Some devices have a <i>Lock-Out</i> mechanism that prevents anyone operating an empty inhaler device, making it clear when a new inhaler is required.</p> <p>The <i>Dose consistency</i> describes the variation in the dose received by the patient, as this often depends on how hard or fast they inhale through the device.</p> <p>The <i>Device Resistance</i> describes how hard it feels when inhaling through this inhaler device. Further explanation of this is described later. This information will be useful to adjust settings of In-Check DIAL inspiratory flow meter (available quarter 2 2016) to measure whether patients can inhale through different inhaler devices at the correct speed of inhalation. This required speed of inhalation to receive a clinically effective dose is described in the <i>Inspiratory Flow</i> section.</p>
Key Features																					
Dexterity	visually required																				
Feedback: Dose Counter																					
Feedback: Tactile/Puff																					
Feedback: Sound																					
Feedback: Visual																					
Dose Lock-Out																					
Dose Consistency																					
Device Resistance	officer resistance																				
Inspiratory Flow																					
<table border="1"> <tr> <td colspan="2">What you need to know / Critical Steps</td> </tr> <tr> <td>Priming</td> <td></td> </tr> <tr> <td>Mouth Position</td> <td></td> </tr> <tr> <td>Dose Warning</td> <td></td> </tr> <tr> <td>Steps</td> <td></td> </tr> <tr> <td>Inspiratory Method</td> <td></td> </tr> <tr> <td>Other Information</td> <td></td> </tr> </table>	What you need to know / Critical Steps		Priming		Mouth Position		Dose Warning		Steps		Inspiratory Method		Other Information		<p>This section describes the essential 'need to know' details to explain the critical steps required to use each device correctly.</p>						
What you need to know / Critical Steps																					
Priming																					
Mouth Position																					
Dose Warning																					
Steps																					
Inspiratory Method																					
Other Information																					

## Types of Inhaler Device

There is a wide variety of different inhalers currently on the market, and can be broadly classified into pressurised metered dose (pMDI), dry powder (DPI), breath-actuated MDI (BA-MDI) and soft mist.

### • Pressurised metered-dose inhalers

pMDIs are widely prescribed as they are cheap and can deliver a wide variety of medications.<sup>4</sup> Despite their wide prescription many patients cannot use pMDIs correctly, even with education and training.<sup>5-7</sup> One study reported that only 79% of patients could use a pMDI correctly after expert training.<sup>8</sup> Many of the difficulties with pMDIs are related to the need to coordinate activation of the device while inhaling slowly and deeply; and approximately 60 and 92% of patients with COPD and asthma, respectively, inhaling too fast from a pMDI.<sup>9,10</sup> The optimal inspiratory flow rate through a pMDI is 25-60 L/min, equivalent to taking a full deep inspiration over 3-5 seconds. When patients are taught to use pMDI devices, they should be advised to inhale **slowly and steadily** through the inhaler.

The use of extra-fine particle pMDIs (e.g. Qvar, Fostair) has been shown to demonstrate equivalent efficacy to other pMDIs but with a reduced equivalent dose of corticosteroid.<sup>11-14</sup> These have been reported to be less dependent upon good inhaler technique than other pMDIs due to longer duration of aerosol emission<sup>15</sup> and lung deposition being affected less by inspiratory flow<sup>16</sup> and coordination.<sup>17</sup>

A spacer is a device that is attached to a pMDI, and produces a reservoir into which the drug aerosol can be generated. This allows the patient to actuate the pMDI without having to coordinate inhalation. There is evidence that the addition of a spacer can increase deposition within the lung of the patient.<sup>18</sup>

### • Breath-actuated metered-dose inhalers

BA-MDIs, such as Autohaler and Easi-Breathe devices, are a development from the original pMDIs. They contain a flow trigger, and the drug is released only when the patient inhales through the inhaler device. BA-MDIs may improve the deposition of aerosol within the lungs of patients who have difficulty coordinating inhalation and actuation of pMDIs.<sup>17,19</sup> When patients are taught to use BA-MDI devices, they should be advised to inhale **slowly and steadily** through the inhaler.

### • Soft mist inhalers

The only current device on the European market is the Respimat. This does not require a propellant or the patient's inspiratory flow to generate an aerosol. The aerosol is generated by a spring that forces the liquid drug formulation through an extremely fine nozzle system, resulting in an aerosol containing a high fine particle fraction.<sup>20</sup> The Respimat has a long aerosol generation time and a low aerosol velocity. These features reduce problems with coordination of actuation and inhalation, and result in a higher lung and lower oropharyngeal deposition.<sup>20</sup> When patients are taught to use a Respimat, they should be advised to inhale **slowly and steadily** through the inhaler.

### • Dry powder inhalers

There are a number of DPIs currently available on the market. Some of these devices are single dose, such as the HandiHaler, Breezhaler or Aerolizer, which require loading of a capsule containing the drug in powder form. Others are multiple dose, such as Turbohaler, Accuhaler, Genuair, Ellipta or NEXThaler.

One of the main advantages that DPIs have over pMDIs are that the generation of the drug aerosol is driven by the patient's inhalation. As a consequence there is no need for an aerosol propellant, or for coordination between actuation of the device and inhalation. DPIs are designed so that force of inhalation creates turbulence within the device, which deaggregates (breaks down) the powder into fine particles that are small enough to be deposited in the lungs.

This turbulence is created by the inspiratory flow generated by the patient through the inhaler device and the internal resistance of the inhaler device itself. This means that low-resistance inhaler DPI devices (such as the Breezhaler) require much faster inspiratory flows than through high-resistance DPI devices (such as the Easyhaler or HandiHaler). If a person cannot generate sufficient inspiratory flow then the drug may not be delivered optimally, as insufficient deaggregation of the dry powder will occur, and large drug particles will be inhaled, which will deposit mainly in the mouth and oropharynx rather than in the lungs. When patients are taught to use any DPI device, they should be advised to inhale as **quickly and deeply** as possible through the inhaler.

### Use of placebo inhalers

For infection control reasons, placebo inhaler devices should be for 'single patient use' only.

Placebo inhaler devices are obtained free-of-charge from manufacturers via their drug reps, or by directly contacting the manufacturer. Ordering information for placebo devices has been collated by the London Medicines Evaluation Network, and is available at: <http://www.medicinesresources.nhs.uk/upload/Availability%20of%20placebo%20inhalers%20FINAL%20June13%20LMEN.pdf>

As it is often difficult to obtain enough placebo inhalers to assess inhaler technique prior to prescribing every new inhaler device, their use is often rationed to use with patients where there is deemed to be a greatest need to teach inhaler technique prior to prescribing inhalers. Decisions about whether or not to use placebo inhalers may be assisted through the use of an In-Check DIAL inspiratory flow meter, which provides an objective measure of a person's to inhale correctly through different devices.

- Placebo devices might be used for patients where it is less certain that they would be able to use an inhaler device. These patients must receive training and assessment *prior* to being given a prescription.
- If patients are judged to be more likely to be able to use an inhaler device, a prescription may be written before inhaler technique training and assessment, which would be performed using their new inhaler devices.

### Use of In-Check DIAL Inspiratory Flow Meters

This is a useful device to aid inhaler technique assessment. It mimics the internal resistance of a range of inhaler devices allowing the measurement of inspiratory flow rate and ensuring the patient can inhale through devices at clinically effective inspiratory flow rates.<sup>21,22</sup>



The current In-Check DIAL 6T&S (Clement Clarke Ltd, Harlow, UK) is validated for just six inhaler devices (Accuhaler, Autohaler, Easi-Breathe, pMDI, Turbohaler [old and new styles]), and so has limited practical value with the wide range of inhaler devices currently available in the UK.

A new, updated version of the In-Check DIAL inspiratory flow meter is anticipated to be launched in the UK in Quarter 2 2016, which will be able to measure inspiratory flow for all inhaler devices, classified by the intrinsic airflow resistance of each device.

The Leeds Inhaler Device Guide documents the airflow resistance of each inhaler device currently available in the UK, and we recommend the use of these new In-Check DIAL inspiratory flow meters for routine inhaler technique assessment, once they are available in the UK.

### Other Tools Available to Assess Inspiratory Flow

#### • Inhaler device whistles



Accuhaler [GlaxoSmithKline]



Ellipta [GlaxoSmithKline]



Turbohaler [AstraZeneca]

Inhaler device whistles are used to assess and ensure patients are able to inhale at, or above, the optimal inspiratory flow necessary to be able to use the device correctly. These are available from the manufacturer.

#### • Flo-tone Trainer [Clement Clarke]

The Flo-Tone trainer is whistle device that attaches to the mouthpiece of a pMDI inhaler device, which is used to teach patients how fast to inhale through pMDI devices. During inhalation, the Flo-Tone will 'whistle', providing an audible signal that the user is inhaling through the pMDI at the correct inspiratory flow and should actuate the inhaler.



### Recommendations for Performing Inhaler Technique Assessments

#### Stepwise Prescribing of Inhaled Medication

1. Confirm the diagnosis
2. Determine the severity of respiratory disease
3. Decide on the class of drug required
4. Teach and assess inhaler technique in order to *Objectively* assess which device the patient can use
  - Is the patient happy with that device?
5. Select the drug that is available in that device

#### Tips for Practical Inhaler Device Selection

1. Inhaler technique should be assessed prior to prescribing inhalers in patients who have never used inhaled medication before.
2. The In-Check DIAL inspiratory flow meter should be used to determine which device(s) the patient is physically able to inhale through.
3. Assess whether your patient is likely to be able to use the device they can achieve the optimal inspiratory flow through. (Take into account age, dexterity, level of understanding, confusion status etc).
  - If likely to be competent, prescribe required drug and train and assess inhaler technique using new inhaler devices.
  - If competency is doubtful or uncertain, train and assess inhaler technique using placebo device.

### Inhaler Technique Assessments

#### Step 1 - Check inspiratory flow

The In-Check DIAL inspiratory flow meter should be used. Depending on the patient and the drug classes required (e.g. SABA, LAMA, LABA, LAMA/LABA, ICS/LABA/etc.), check the patient's inspiratory flow through a range of different devices. It is useful to check inspiratory flow two or three times through each device to ensure that the patient achieves the same inspiratory flows consistently.

#### Step 2 - Teach correct inhaler technique

Where possible, education on correct inhaler technique requires the healthcare professional to demonstrate correct inhaler technique to the patient. It will be useful to keep your own set of placebos for you to use.

- Data from a small study in asthmatic patients identified that the most successful method of teaching correct inhaler technique is to demonstrate correct inhaler technique using placebo inhalers so that the patient can copy the demonstrator, and to explain each step with emphasis on the critical steps (those that if the patient fails to perform they receive no dose e.g. failing to remove the cap on a pMDI).<sup>23</sup>

#### Step 3 - Check inhaler technique and check understanding

Once the patient has been taught and shown how to use their inhaler, they should be asked to demonstrate how they would use it. This allows the healthcare professional to check they understand how to use their inhaler device, and to reinforce any steps they are unable to perform correctly.

#### Step 4 - Educate patient on rationale for using inhaled medication

Adherence to prescribed inhaled medication may be improved if the rationale for using each medicine is explained to the patient. This should also include when and how to use their inhaler devices.

## Aerosol science & inhaled medication

### • Effects of particle size on lung deposition

Depending on their particle size, inhaled drug particles will deposit in different regions of the lung. Particles  $<1 \mu\text{m}$  are likely to reach the peripheral airways and alveoli or will be exhaled, particles  $1\text{--}5 \mu\text{m}$  will deposit in the large and conducting airways, while particles  $>5 \mu\text{m}$  will predominately deposit in the oropharynx.<sup>24,25</sup>

There are three different mechanisms of deposition of aerosols: inertial impaction, sedimentation and diffusion (figure 1). However for the particle size used in aerosol therapies of approximately  $1\text{--}10 \mu\text{m}$ , only two of these mechanisms predominate: inertial impaction and gravitational sedimentation. The third mechanism, Brownian motion/diffusion is only relevant in aerosols less than 1 micron in diameter, and consequently is unlikely to be important for inhaled drugs.

Inertial impaction occurs in either the oropharynx, or at bifurcations of main branches of the bronchial tree, particularly in the large central airways. It occurs mainly with large particles or high velocity particles (i.e. those with high inertia), where they are unable to follow the airstream when it changes direction, thus impacting on the airway wall.<sup>24,26</sup>

Gravitational sedimentation occurs for smaller particles that are able to follow the airstream and penetrate the more peripheral bronchioles and alveoli. Here, the airstream flows slower, allowing the particles to settle on to the airway surfaces either during the course of slow steady breathing or during breath-holding.<sup>24,26</sup> Breath-holding is important for smaller particle sizes owing to the increased chance of exhalation of the drug, because they can remain airborne for a considerable time.<sup>16</sup>




Figure 1. Particle deposition in the respiratory tract. (Adapted with permission)<sup>26</sup>

### Implications for practice:

Particle size is important: those that are too small may be exhaled; those that are too large experience inertial impaction in the oropharynx and large conducting airways.

Increased aerosol particle speed increases the probability of deposition by impaction in the oropharynx and large conducting airways; slow aerosol particle speed allows more particles to penetrate the peripheral bronchial tree.

Increasing the inhalation volume allows the aerosol to penetrate peripheral bronchioles.

Breath-holding increases gravitational sedimentation.

## Effects of inspiratory flow on lung deposition

The total lung deposition of an inhaled drug is strongly affected by the speed of inhalation. DPIs require a fast and deep inhalation to 'suck up' the drug in the inhaler device. A fast inhalation rate generates a large internal turbulent force in the inhaler device, which is required to break up the formulation of the metered dose to produce particles of a size distribution that will penetrate the peripheral airways.<sup>26,27,28</sup> Failure to achieve this high internal force increases the likelihood of the dose impacting in the mouth and throat. By contrast, aerosol inhalers, such as the MDI, require a slow and deep inhalation, with an inspiratory flow rate of less than 60 l/min. This is owing to the device generating its own aerosol, and so a slower inhalation rate is required to ensure that the drug deposits in the peripheral airways, since a fast inhalation will increase the velocity of the drug particles, thus increasing inertial impaction in the oropharynx as described above.<sup>10,27</sup> Some of the effect of excessive inspiratory flow with a pMDI can be mitigated by using pMDI devices with smaller particle sizes in the aerosol, which may have a greater lung deposition if the inspiratory flow from 30.8 to 67.1 L/min.<sup>16</sup> However these faster inspiratory flows will reduce lung deposition of  $6\text{-}\mu\text{m}$  particles.<sup>16</sup>

Faster inspiratory flow rates through a pMDI (180 vs. 30 L/min) increases inertial impaction of aerosols in the oropharynx and at bifurcations in the large central airways, thus reducing lung deposition in the peripheral airways by 33%.<sup>29</sup> A significant proportion of patients with asthma and COPD have been shown to have an inspiratory flow that is much too high for an MDI, which may reduce the clinical effectiveness of inhaled drugs.<sup>9,10</sup>

Other studies have demonstrated that lung deposition from DPI devices increases as inspiratory flow increases. This is particularly so for Turbohaler, and to a lesser degree for Accuhaler, devices where lung deposition increases with faster inspiratory flows.<sup>29,30</sup>

As a consequence, when a patient uses an aerosol inhaler such as the MDI, a slow inspiratory flow rate will produce significantly greater lung deposition than when used at a faster inspiratory flow rate. In contrast, when dry powder inhalers such as a Turbohaler or Accuhaler are used, a faster inspiratory flow rate will produce significantly greater lung deposition than a slower inspiratory flow rate.

The ability of patients to use specific inhalers can vary with time. During times of increased hyperinflation, such as acute exacerbations of asthma, or progression to severe emphysema, patients may not be able to generate sufficient inspiratory flow through high-resistance devices. In this situation, a low-resistance device, such as a pMDI or pMDI plus spacer, may be more appropriate.


## References

1. National Institute for Health and Clinical Excellence. NICE Guideline CG101: Chronic obstructive pulmonary disease in over 16s: diagnosis and management. 2010. Available at: <http://www.nice.org.uk/guidance/CG101>
2. Global Initiative for Chronic Obstructive Lung Disease. Global strategy for the diagnosis, management and prevention of COPD. 2015. Available at [www.goldcopd.com](http://www.goldcopd.com).
3. British Thoracic Society/Scottish Intercollegiate Guidelines Network. British Guideline on the management of asthma. 2014. Available at: <https://www.brit-thoracic.org.uk/document-library/clinical-information/asthma/btssign-asthma-guideline-2014/>
4. Lavorini F, Corrigan CJ, Barnes PJ *et al*. Retail sales of inhalation devices in European countries: so much for a global policy. *Respir Med* 2011;105:1099–103.
5. Crompton GK. Problems patients have using pressurized aerosol inhalers. *Eur J Respir Dis* 1982;119 (Suppl.):101–4.
6. Allen SC *et al* Prior A. What determines whether an elderly patient can use a metered dose inhaler correctly? *Br J Dis Chest* 1986;80:45–9.
7. Chapman KR, Love L *et al* Brubaker H. A comparison of breath-actuated and conventional metered-dose inhaler inhalation techniques in elderly subjects. *Chest* 1993;104:1332–7.
8. Lenney J, Innes JA *et al* Crompton GK. Inappropriate inhaler use: assessment of use and patient preference of seven inhalation devices. *Respir Med* 2000;94:496–500.
9. Al-Showair RA, Tarsin WY, Assi KH *et al*. Can all patients with COPD use the correct inhalation flow with all inhalers and does training help? *Respir Med* 2007;101:2395–401.
10. Al-Showair RA, Pearson SB *et al* Chrystyn H. The potential of a 2Tone Trainer to help patients use their metered-dose inhalers. *Chest* 2007;131:1776–82.
11. Fireman P, Prenner BM, Vincken W *et al*. Long-term safety and efficacy of a chlorofluorocarbon-free beclomethasone dipropionate extrafine aerosol. *Ann Allergy Asthma Immunol* 2001;86:557–65.
12. Gross G, Thompson PJ, Chervinsky P *et al*. Hydrofluoroalkane-134a beclomethasone dipropionate, 400 microg, is as effective as chlorofluorocarbon beclomethasone dipropionate, 800 microg, for the treatment of moderate asthma. *Chest* 1999;115:343–51.
13. Papi A, Paggiaro P, Nicolini G *et al*. Beclomethasone/formoterol vs fluticasone/salmeterol inhaled combination in moderate to severe asthma. *Allergy* 2007;62:1182–8.
14. Papi A, Paggiaro PL, Nicolini G *et al*. Beclomethasone/formoterol versus budesonide/formoterol combination therapy in asthma. *Eur Respir J* 2007;29:682–9.
15. Fabbri LM, Nicolini G, Olivieri D *et al*. Inhaled beclomethasone dipropionate/formoterol extra-fine fixed combination in the treatment of asthma: evidence and future perspectives. *Expert Opin Pharmacother* 2008;9:479–90.
16. Usmani OS, Biddiscombe MF *et al* Barnes PJ. Regional lung deposition and bronchodilator response as a function of beta2-agonist particle size. *Am J Respir Crit Care Med* 2005;172:1497–1504.
17. Leach CL, Davidson PJ, Hasselquist BE *et al*. Influence of particle size and patient dosing technique on lung deposition of HFA-beclomethasone from a metered dose inhaler. *J Aerosol Med* 2005;18:379–85.
18. Newman SP. Spacer devices for metered dose inhalers. *Clin Pharmacokinet* 2004;43:349–60.
19. Newman SP, Weisz AW, Talaei N *et al*. Improvement of drug delivery with a breath actuated pressurised aerosol for patients with poor inhaler technique. *Thorax* 1991;46:712–6.
20. Dalby R, Spallek M *et al* Voshaar T. A review of the development of Respimat Soft Mist Inhaler. *Int J Pharm* 2004;283:1–9.
21. Nsour WM, Alldred A, Corrado OJ *et al*. Measurement of peak inhalation rates with an In-Check® Meter to identify an elderly patient's ability to use a Turbuhaler®. *Respir Med* 2001;95:965–8.
22. Chrystyn H. Is inhalation rate important for a dry powder inhaler? Using the In-check Dial to identify these rates. *Respir Med* 2003;97:181–7
23. Basheti IA, Reddel HK *et al* Armour CL. Counseling about turbuhaler technique: needs assessment and effective strategies for community pharmacists. *Respir Care* 2005;50:617–23.
24. Newman SP. Aerosol deposition considerations in inhalation therapy. *Chest* 1985;88(Suppl. 2):152S–160S.
25. Everard ML. Guidelines for devices and choices. *J Aerosol Med* 2001;14(Suppl. 1):S59–S64.
26. Bell J. Why optimise inhaler technique in asthma and COPD? *Br J Prim Care Nurs* 2008;2:37–9.
27. Haughney J, Price D, Barnes NC *et al*. Choosing inhaler devices for people with asthma: current knowledge and outstanding research needs. *Respir Med* 2010;104:1237–45.
28. Chrystyn H. Is inhalation rate important for a dry powder inhaler? Using the In-Check Dial to identify these rates. *Respir Med* 2003;97:181–7.
29. Borgstrom L, Bondesson E, Moren F *et al*. Lung deposition of budesonide inhaled via Turbuhaler: a comparison with terbutaline sulphate in normal subjects. *Eur Respir J* 1994;7:69–73.
30. Palander A, Mattila T, Karhu M *et al*. In vitro comparison of three salbutamol containing multidose dry powder inhalers: Buventol Easyhaler®, Inspieryl Turbuhaler® and Ventoline Diskus®. *Clin Drug Investig* 2000;20:25–33.

### This section is adapted from:

- Capstick TGD *et al* Clifton IJ. Inhaler technique and training in people with chronic obstructive pulmonary disease and asthma. *Expert Rev Respir Med* 2012;6:91–103.
- Dimov D. COPD - Full Guidelines for the Diagnosis and Management of Chronic Obstructive Pulmonary Disease. Last Review: 18/03/2013. Available at: <http://www.lhp.leedsth.nhs.uk/common/guidelines/detail.aspx?ID=572>

## ACCUHALER




<b>Type of inhaler:</b>	Single dose DPI - doses contained in separate blisters in a foil strip
<b>In use expiry:</b>	No restrictions
<b>Available as</b>	
<b>SABA</b>	Ventolin (salbutamol)
<b>SAMA</b>	
<b>LABA</b>	Serevent (salmeterol)
<b>LAMA</b>	
<b>LAMA/LABA</b>	
<b>ICS</b>	Flixotide (fluticasone)
<b>ICS/LABA</b>	Seretide (fluticasone/salmeterol)

Photo. ©GlaxoSmithKline

Key Features	
<b>Dexterity</b>	Low-Moderate dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts individual doses. Last 5 numbers highlighted in red text
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	No lock-out - lever can continue to be pushed back when device is empty
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (30-90 L/min)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	Inspiratory flow - test with In Check DIAL inspiratory flow meter Inspiratory flow - test with Accuhaler whistle

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use, and then open inhaler.
<b>2. Priming</b>	Hold inhaler horizontally, or with mouthpiece uppermost. <b>Push lever back completely</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cover after use.


What you need to know / Critical Steps	
<b>Priming</b>	Ensure lever is pushed back completely. Do not close device before inhalation
<b>Multi-Dosing</b>	Risk of multi-dosing from multiple actuations
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Inhaler may sound/feel gritty if incomplete inhalation of dose
<b>Cleaning</b>	Wipe the mouthpiece with a dry tissue

AEROLISER	
 <p>Photo. ©Novartis</p>	<b>Type of inhaler:</b> Single dose DPI - doses contained in separate capsules
	<b>In use expiry:</b> No restrictions
	<b>Available as</b>
	<b>SABA</b>
	<b>SAMA</b>
	<b>LABA</b> Foradil (formoterol)
	<b>LAMA</b>
<b>LAMA/LABA</b>	
<b>ICS</b>	
<b>ICS/LABA</b>	

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	No dose counter. Placement of capsule for each dose.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	Whirring sound during inhalation
<b>Feedback: Visual</b>	Transparent capsule should be empty
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Remove the cap and rotate the mouthpiece of the inhaler towards the arrow on the bottom of the inhaler to open.
<b>2. Priming</b>	Remove one capsule from the blister and place into the chamber in the middle of the inhaler. Twist the mouthpiece in the opposite direction to close, until you hear a click. Press the side buttons inward to pierce the capsule and release.
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply</b> - you should hear a whirring noise
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Remove the capsule and replace cap.


What you need to know / Critical Steps	
<b>Priming</b>	Capsule must be placed into the chamber and pierced.
<b>Multi-Dosing</b>	Single dose - use one capsule at a time.
<b>Dose Wasting</b>	Wasted dose if inhalation not repeated and capsule not emptied.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Strong and deep
<b>Other information</b>	Occasionally small fragments of the capsule may be inhaled (usually if the capsule is pierced more than once). This is not harmful. Capsules can get stuck in the inhaler and no 'whirring' noise will be heard during inhalation. If this occurs tap the base of inhaler on a hard surface to release the capsule.
<b>Cleaning</b>	Wipe the mouthpiece and capsule compartment with a dry cloth or a soft clean brush .

AUTOHALER		
 <p>Photo. ©Teva</p>	<b>Type of inhaler:</b> Breath-actuated MDI	
	<b>In use expiry:</b> No restrictions	
	<b>Available as</b>	
	SABA	Airimir (salbutamol)
	SAMA	
	LABA	
	LAMA	
LAMA/LABA		
ICS	Qvar (beclometasone)	
ICS/LABA		

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	No dose counter
<b>Feedback: Taste/Feel</b>	Taste of aerosol. Both Autohalers contain small amounts of ethanol (maybe unpleasant taste for children)
<b>Feedback: Sound</b>	Click on inhalation
<b>Feedback: Visual</b>	Lever should be in upright position before inhalation and lowered down after each puff.
<b>Device Lock-Out</b>	No lock-out - lever can continue to be pushed up when device is empty
<b>Dose Consistency</b>	Consistent metered dose across a range of inspiratory flow rates (30-60 L/min); significantly reduced at fast inspiratory flow rates (>60L/min)
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	Test with In Check DIAL inspiratory flow meter

Instructions for Use	
<b>1. Preparation</b>	Remove the cap
<b>2. Priming</b>	Shake the inhaler. Keeping the inhaler upright, without blocking
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale slowly and deeply - do not stop when device 'clicks'
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 2-6 to ensure you inhaled the full dose. Replace cap after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position. Ensure the lever is pushed upright.
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Wasted dose if dose release slide used at the bottom of the inhaler
<b>Vents</b>	Air vents at the bottom of the device
<b>Inspiratory Method</b>	Slow and steady
<b>Other information</b>	Lever must be pushed into upright position before each dose.
<b>Cleaning</b>	Clean the mouthpiece weekly with a dry tissue or cloth

BREEZHALER	
 <p>Photo. ©Novartis</p>	<b>Type of inhaler:</b> Single dose DPI - doses contained in separate capsules
	<b>In use expiry:</b> 30 days after first use
	<b>Available as</b>
	<b>SABA</b>
	<b>SAMA</b>
	<b>LABA</b> Onbrez (indacaterol)
	<b>LAMA</b> Seebri (glycopyrronium)
<b>LAMA/LABA</b> Ultibro (glycopyrronium/indacaterol)	
<b>ICS</b>	
<b>ICS/LABA</b>	

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	No dose counter. Placement of capsule for each dose.
<b>Feedback: Taste/Feel</b>	Taste of powder
<b>Feedback: Sound</b>	Whirring sound during inhalation
<b>Feedback: Visual</b>	Transparent capsule should be empty
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (50-100 L/min)
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Remove the cap and tilt the mouthpiece to open the inhaler.
<b>2. Priming</b>	Remove one capsule from the blister and place into the chamber in the middle of the inhaler. Close the mouthpiece and press down until you hear a click. Press the side buttons inward to pierce the capsule and release.
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale strongly and deeply - you should hear a whirring noise
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Remove the capsule and replace cap.

What you need to know / Critical Steps	
<b>Priming</b>	Capsule must be placed into the chamber and pierced.
<b>Multi-Dosing</b>	Single dose - use one capsule at a time.
<b>Dose Wasting</b>	Wasted dose if inhalation not repeated and capsule not emptied.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Strong and deep
<b>Other information</b>	Occasionally small fragments of the capsule may be inhaled (usually if the capsule is pierced more than once). This is not harmful. Capsules can get stuck in the inhaler and no 'whirring' noise will be heard during inhalation. If this occurs tap the base of inhaler on a hard surface to release the capsule.
<b>Cleaning</b>	Clean the mouthpiece inside and outside weekly with a dry tissue or lint-free cloth


CLICKHALER	
	<b>Type of inhaler:</b> Reservoir Multidose DPI
	<b>In use expiry:</b> 6 months after removed from foil pouch
<b>Available as</b>	
SABA	Asmasal (salbutamol)
SAMA	
LABA	
LAMA	
LAMA/LABA	
ICS	Asmabec (beclometasone)
ICS/LABA	

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Dose counter
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	Device locks when empty
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (15-60 L/min)
<b>Device Resistance</b>	Medium airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Remove the cap.
<b>2. Priming</b>	Shake the inhaler, and then hold it upright, with your thumb on the base. <b>Press the button down firmly, until you hear a click, then release.</b>
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cover after use.

What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the upright position
<b>Multi-Dosing</b>	Risk of multi-dosing from multiple actuations
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting or inverting the device may lose the dose. Wasted dose if multiple actuations used during a single breath
<b>Vents</b>	Air vents in the middle of the device
<b>Inspiratory Method</b>	Strong and deep
<b>Other information</b>	No coordination is required.
<b>Cleaning</b>	Press down on the end of the mouthpiece to remove it and wipe with a dry cloth or tissue.

## EASI-BREATHE




<b>Type of inhaler:</b>	Breath-actuated MDI
<b>In use expiry:</b>	No restrictions
<b>Available as</b>	
SABA	Salamol (salbutamol)
SAMA	
LABA	
LAMA	
LAMA/LABA	
ICS	Qvar (beclometasone)
ICS/LABA	

Photo. ©Teva

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	No dose counter
<b>Feedback: Taste/Feel</b>	Taste of aerosol
<b>Feedback: Sound</b>	Audible aerosol spray when dose delivered
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent metered dose across a range of inspiratory flow rates (20-60 L/min); significantly reduced at fast inspiratory flow rates (>60L/min)
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	Test with In Check DIAL inspiratory flow meter

Instructions for Use	
<b>1. Preparation</b>	Remove the cap
<b>2. Priming</b>	Shake the inhaler, and then hold it upright, with your thumb on the base.
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale slowly and deeply- do not stop when device 'puffs'
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Replace cap after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the upright position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	n/a
<b>Vents</b>	Do not block air vents at the top of the device
<b>Inspiratory Method</b>	Slow and deep
<b>Other information</b>	Do not unscrew top half of Easi-Breathe device
<b>Cleaning</b>	Salamol: Unscrew and remove the top of the inhaler and remove the metal cannister. Rinse the bottom of the empty inhaler in warm running water for 30 seconds. Then dry thoroughly (leave overnight if possible) and reassemble.  Qvar: Clean the mouthpiece weekly with a dry tissue or cloth.

EASYHALER															
	<p><b>Type of inhaler:</b> Reservoir Multidose DPI</p> <p><b>In use expiry:</b> 6 (salbutamol, beclometasone and budesonide) and 4 (formoterol) months after removing from foil pouch</p> <p><b>Available as</b></p> <table border="1"> <tr> <td>SABA</td> <td>Easyhaler salbutamol</td> </tr> <tr> <td>SAMA</td> <td></td> </tr> <tr> <td>LABA</td> <td>Easyhaler formoterol</td> </tr> <tr> <td>LAMA</td> <td></td> </tr> <tr> <td>LAMA/LABA</td> <td></td> </tr> <tr> <td>ICS</td> <td>Easyhaler beclometasone; Easyhaler budesonide</td> </tr> <tr> <td>ICS/LABA</td> <td></td> </tr> </table>	SABA	Easyhaler salbutamol	SAMA		LABA	Easyhaler formoterol	LAMA		LAMA/LABA		ICS	Easyhaler beclometasone; Easyhaler budesonide	ICS/LABA	
SABA	Easyhaler salbutamol														
SAMA															
LABA	Easyhaler formoterol														
LAMA															
LAMA/LABA															
ICS	Easyhaler beclometasone; Easyhaler budesonide														
ICS/LABA															
Photo. ©Orion															

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts down in steps of 10 doses. Red indicator begins to show with 20 doses remaining.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent dosing across range of inspiratory flow rates (30-60 L/min)
<b>Device Resistance</b>	High airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use, and then remove cap.
<b>2. Priming</b>	Shake the inhaler, then <b>whilst holding the inhaler upright press the button all the way down and then release it.</b> You should hear a click sound.
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Unused doses collect in bottom chamber with clear window
<b>Cleaning</b>	Clean the mouthpiece weekly with a dry tissue or cloth.

ELLIPTA	
 <p>Photo. ©GlaxoSmithKline</p>	<b>Type of inhaler:</b> Single dose DPI - doses contained in separate blisters in a foil strip
	<b>In use expiry:</b> 6 weeks after opening the tray
	<b>Available as</b>
	<b>SABA</b>
	<b>SAMA</b>
	<b>LABA</b>
	<b>LAMA</b>
<b>LAMA/LABA</b>	Increase (umeclidinium)
<b>ICS</b>	Anoro (umeclidinium/vilanterol)
<b>ICS/LABA</b>	Relvar (fluticasone furoate/vilanterol)

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts down each dose. Red indicator when less than 10 doses remain.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent dosing across range of inspiratory flow rates (43.5-130 L/min)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	No validated test available. Test with Ellipta whistle.

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use.
<b>2. Priming</b>	Hold inhaler in the upright position. <b>Open the cover completely.</b>
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale strongly and deeply.
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cover after use.

What you need to know / Critical Steps	
<b>Priming</b>	Ensure cover is pushed back completely. Do not close device before inhalation
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Potential to waste doses if cap is opened and closed more than once without inhalation. Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	Do not block air vents
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	No coordination is required.
<b>Cleaning</b>	Use a dry tissue to clean the mouthpiece

FORSPIRO	
 <p>Photo. ©Sandoz</p>	<b>Type of inhaler:</b> Single dose DPI - doses contained in separate blisters in a foil strip
	<b>In use expiry:</b> No restrictions
	<b>Available as</b>
	<b>SABA</b>
	<b>SAMA</b>
	<b>LABA</b>
	<b>LAMA</b>
<b>LAMA/LABA</b>	
<b>ICS</b>	
<b>ICS/LABA</b> AirFluSal (fluticasone/salmeterol)	

Key Features	
<b>Dexterity</b>	Moderate-high dexterity required. Requires significant manipulation.
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts individual doses. Last 10 doses highlighted with red border.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	Loading and priming of each dose can be observed. Used foil blisters will appear in the side chamber (although used blisters will not appear until 2 days after each one is used)
<b>Device Lock-Out</b>	No lock-out - lever can continue to be pushed back when device is empty
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (30-90 L/min)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Ensure side chamber is closed. <b>Open cap, then check dose counter before use.</b>
<b>2. Priming</b>	Hold inhaler horizontally, or with mouthpiece uppermost. <b>Push white lever open completely until it clicks, then close it again fully until it clicks back into original position</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for 5 to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cap after use.

What you need to know / Critical Steps	
<b>Priming</b>	Ensure white lever is pushed open completely until it clicks, then closed again fully until it clicks back into original position.
<b>Multi-Dosing</b>	Single dose - reportedly cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	Do not block air inlets to the side of the mouthpiece
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	The foil strip must be torn away against the 'teeth' of the side chamber. Failure to do so may cause the inhaler to jam.
<b>Cleaning</b>	Clean the outside of the mouthpiece with a dry tissue



GENUAIR	
	<b>Type of inhaler:</b> Reservoir Multidose DPI
	<b>In use expiry:</b> 90 (Eklira) or 60 (Duaklir) days after removal from foil pouch
	<b>Available as</b>
	SABA
	SAMA
	LABA
	LAMA Eklira (aclidinium)
	LAMA/LABA Duaklir (aclidinium/formoterol)
ICS	
ICS/LABA	

Photo. ©AstraZeneca

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Dose counter. Counts down in steps of 10 doses. Red indicator shows for last 10 doses.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	'Click' on inhalation, and 'whirring' noise during inhalation
<b>Feedback: Visual</b>	Coloured control window changes from red to green after priming (pressing the button), then back to red after inhalation
<b>Device Lock-Out</b>	Device locks when empty
<b>Dose Consistency</b>	Consistent dosing across range of inspiratory flow rates (so long as inhalation speed >35L/min to cause device to click)
<b>Device Resistance</b>	Medium airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Check dose counter and ensure that the coloured control window is red before use. Remove the cap.
<b>2. Priming</b>	Hold the inhaler upright with the button facing up. Press the button all the way down and then release it.
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale strongly and deeply - do not stop when device 'clicks'
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Device will not 'click' with low inhalation speed (low emitted dose), preventing further doses to be primed
<b>Cleaning</b>	Wipe the mouthpiece with a dry tissue.

HANDIHALER	
	<b>Type of inhaler:</b> Single dose DPI - each dose contained in separate capsules.
	<b>In use expiry:</b> Capsules: 9 days after opening the blister. HandiHaler device: discard after 12 months
	<b>Available as</b>
	SABA
	SAMA
	LABA
	LAMA
LAMA/LABA	
ICS	
ICS/LABA	
Spiriva (tiotropium)	

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	Remaining capsules may be counted.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose). Capsule will vibrate during inhalation
<b>Feedback: Sound</b>	Vibration of capsule heard and/or felt during inhalation
<b>Feedback: Visual</b>	Open to visualise empty capsule after use.
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent dosing over low to moderate inspiratory flow (28-60 L/min), but significant reduction at lower inspiratory flows.
<b>Device Resistance</b>	High airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Open dust cap and mouthpiece
<b>2. Priming</b>	Remove capsule from blister strip and insert capsule into centre chamber. Close mouthpiece until a click is heard and keeping the inhaler upright, press the piercing button once and release.
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale slowly and deeply, but at a rate sufficient to hear the capsule vibrate.
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Open mouthpiece to empty the used capsule and close the mouthpiece and dust cap.


What you need to know / Critical Steps	
<b>Priming</b>	Capsule must be placed into the chamber and pierced.
<b>Multi-Dosing</b>	Single dose - cannot multi-dose.
<b>Dose Wasting</b>	Exhaling into device will lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Occasionally small fragments of the capsule may be inhaled (usually if the capsule is pierced more than once). This is not harmful.
<b>Cleaning</b>	Clean monthly. Open the dust cap and mouthpiece, and also the base by lifting the piercing button. Rinse the inhaler with warm water to remove any powder. Leave to dry for 24 hours. If needed, the mouthpiece may be cleaned with a moist tissue.

NEXTHALER		
 <p>Photo. ©Chiesi</p>	<b>Type of inhaler:</b> Reservoir Multidose DPI	
	<b>In use expiry:</b> 6 months after removed from foil pouch	
	<b>Available as</b>	
	SABA	
	SAMA	
	LABA	
	LAMA	
LAMA/LABA		
ICS		
ICS/LABA	Fostair (beclometasone/formoterol)	

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts down each dose. Red arrows indicate when less than 10 doses remain.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	'Click' on inhalation
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (30-90 L/min)
<b>Device Resistance</b>	Medium-high airflow resistance
<b>Inspiratory Flow</b>	No validated test available. Breath actuated mechanism -device will only release dose with correct speed of inhalation.

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use.
<b>2. Priming</b>	Hold inhaler in the upright position, and shake the device. Open the cover completely.
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale strongly and deeply - do not stop when device 'clicks'.
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cover after use.

What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position.
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Contains dose protector - exhaling into device will NOT lose the dose
<b>Vents</b>	Do not block air vents
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Device will not 'click' with low inhalation speed (no dose will be released).
<b>Cleaning</b>	Use a dry tissue or cloth to clean the inhaler

NOVOLIZER	
	<b>Type of inhaler:</b> Reservoir Multidose DPI
	<b>In use expiry:</b> Cartridge: 3 months after opening. Novolizer device: 1 year or 2,000 doses (which ever comes first)
	<b>Available as</b>
SABA	Salbutin (salbutamol)
SAMA	
LABA	
LAMA	
LAMA/LABA	
ICS	Budelin (budesonide)
ICS/LABA	

### Instructions for Use


**Inserting the Cartridge**  
 1. Lightly press together the ribbed buttons on the lid, and pull gently forwards to lift it off.  
 2. Take the new cartridge out of the aluminium foil packaging.  
 3. Insert the cartridge into the Novolizer with the dose counter facing the mouthpiece.  
 4. Replace the lid on the Novolizer device.

<b>1. Preparation</b>	<b>Check dose counter and ensure that the coloured control window is red before use. Remove the cap.</b>
<b>2. Priming</b>	<b>Hold the inhaler upright with the button facing up. Press the button all the way down and then release it.</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply</b> - do not stop when device 'clicks'
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use.

Key Features	
<b>Dexterity</b>	Low-moderate dexterity required
<b>Feedback: Dose Counter</b>	Dose counter. Counts down in steps of 20 doses. Shaded indicator shows for last 20 doses.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	'Click' on inhalation, and 'whirring' noise during inhalation
<b>Feedback: Visual</b>	Coloured control window changes from red to green after priming (pressing the button), then back to red after inhalation
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Inconsistent dosing across range of inspiratory flow rates (lower effective dose with low inhalation speed)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

### What you need to know / Critical Steps

<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Device will not 'click' with low inhalation speed (low emitted dose), preventing further doses to be primed
<b>Cleaning</b>	Clean regularly and when the cartridge is changed. Remove the cap and then the mouthpiece by twisting it anticlockwise. Turn the device upside down and remove the lid. Tap to remove any remaining powder. Clean the mouthpiece, dispensing slide and powder inhaler with a soft, lint-free dry cloth.

pMDI		
	<b>Type of inhaler:</b> Pressurised Metered Dose Inhaler	
	<b>In use expiry:</b> No restrictions, except: Atimos Modulite: 3 months at room temp.; Flutiform: 3 months after removal from foil pouch; Fostair: 5 months at room temp.	
	<b>Available as</b>	
	<b>SABA</b>	Airomir, AirSalb, Salamol, Ventolin (salbutamol)
	<b>SAMA</b>	Atrovent (ipratropium)
	<b>LABA</b>	Atimos Modulite (formoterol); Neovent, Serevent, Vertine (salmeterol)
	<b>LAMA</b>	
	<b>LAMA/LABA</b>	
<b>ICS</b>	Clenil Modulite, Qvar (beclometasone); Flixotide (fluticasone)	
<b>ICS/LABA</b>	Flutiform (fluticasone/formoterol); Fostair (beclometasone/formoterol); Seretide, Sirdupla (fluticasone/salmeterol)	

Instructions for Use	
<b>1. Preparation</b>	Remove the cap.
<b>2. Priming</b>	Shake the MDI inhaler, and then hold it upright, with your thumb on the base.
<b>3. Exhaling</b>	Exhale fully and away from the mouthpiece.
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale slowly and deeply, pressing the canister down to release the medicine at the start of inhalation, and continue to inhale deeply</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Replace cap after use.

Key Features	
<b>Dexterity</b>	Moderate dexterity required. Haleraid available for GSK MDIs to assist actuation for people with physical limitations
<b>Feedback: Dose Counter</b>	Usually no dose counter. (Exceptions: Flutiform, Seretide, Sirdupla - counts each dose)
<b>Feedback: Taste/Feel</b>	Taste of aerosol. Some MDIs contain small amounts of ethanol (maybe unpleasant taste for children)
<b>Feedback: Sound</b>	Audible aerosol spray
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent metered dose across a range of inspiratory flow rates (25-60 L/min); significantly reduced at fast inspiratory flow rates (>60L/min)
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	Inspiratory flow - test with In Check DIAL inspiratory flow meter Inspiratory flow - test with Flo-Tone whistle & Trainhaler

What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Wasted dose if multiple actuations used during a single breath
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Slow and steady
<b>Other information</b>	Good co-ordination of inhalation and actuation is required
<b>Cleaning</b>	Clean the mouthpiece with a dry tissue or cloth

## pMDI + SPACER (Multiple breath method)

Available Spacers  
Volumatic ± mask



<b>Type of inhaler:</b>	Pressurised Metered Dose Inhaler + large volume spacer
<b>In use expiry:</b>	Replace after 6 to 12 months
<b>Compatible MDIs</b>	Ensure MDI is compatible with spacer before prescribing. Volumatic Spacer is generally only compatible with GSK MDIs
<b>When to use</b>	Useful for: <ul style="list-style-type: none"> <li>• Infants and toddlers.</li> <li>• Reducing side effects with high doses of inhaled corticosteroids.</li> <li>• People with poor co-ordination.</li> <li>• People who are unable to hold their breath after inhalation, and during exacerbations.</li> </ul>

### Key Features


<b>Dexterity</b>	Low-moderate dexterity required
<b>Feedback: Dose Counter</b>	n/a
<b>Feedback: Taste/Feel</b>	May reduce taste of pMDI aerosol
<b>Feedback: Sound</b>	Click of valve audible when breathing in and out
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent dosing when used correctly.
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

### Instructions for Use

<b>1. Preparation</b>	Remove the cap from inhaler, and assemble the two parts of the Volumatic spacer.
<b>2. Priming</b>	Shake the inhaler, then <b>place the mouthpiece of the inhaler into flat end of the spacer</b> and hold it upright, with your thumb on the base
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	<b>Place mouthpiece of spacer between teeth and lips and actuate one dose into spacer</b>
<b>5. Inhalation</b>	<b>Inhale slowly and deeply through the mouth for 5 to 6 breaths.</b>
<b>6. Breath holding</b>	Then remove the mouthpiece of the spacer from your mouth. No breath-holding manoeuvre is required.
<b>7. Closing &amp; Repeating</b>	Remove inhaler from the spacer and replace cap.

### What you need to know / Critical Steps


<b>Priming</b>	Delay between actuation of MDI into Spacer and inhalation results in reduced effective dose
<b>Multi-Dosing</b>	Multiple actuations of MDI into Spacer results in large particle sizes and reduced effective dose.
<b>Dose Wasting</b>	Important only to use single actuations into the spacer.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Slow and steady. Multiple breath
<b>Other information</b>	Useful for people with poor co-ordination. Useful when high doses of inhaled corticosteroid are prescribed
<b>Cleaning</b>	Clean twice a week. Take the two halves apart and wash with warm water containing a mild detergent, using a soft toothbrush or bottle brush if necessary. Do not rub with a cloth as this creates static. Rinse with clean water and dry at room temperature. Ensure the valve runs freely prior to re-use.

pMDI + SPACER (Single breath method)			
<b>Available Spacers:</b> <ul style="list-style-type: none"> <li>• AeroChamber Plus (+/- mask)</li> <li>• AeroChamber Plus for Infant or Toddler (+ mask)</li> <li>• Able Spacer (+/- mask)</li> <li>• A2A (+/- mask)</li> <li>• DispozABLE Spacer</li> <li>• Optichamber Advantage</li> <li>• Optichamber Diamond (+ mask)</li> <li>• Pocket Chamber (+/- mask)</li> <li>• Space Chamber Plus</li> <li>• Volumatic Spacer (+/- mask)</li> <li>• Vortex (+/- mask)</li> </ul> 	<b>Type of inhaler:</b> Pressurised Metered Dose Inhaler + small volume spacer		
	<b>In use expiry:</b>	Check manufacturer's instructions. AeroChamber Plus may need to be replaced after 12 months. DispozABLE Spacer - single use for emergencies, such as ambulances and schools. Volumatic spacer should be replaced after 6 to 12 months.	
	<b>Compatible MDIs</b>	Ensure MDI is compatible with spacer before prescribing.	
	<b>When to use</b>	Useful for: <ul style="list-style-type: none"> <li>• Infants and toddlers.</li> <li>• Reducing side effects with high doses of inhaled corticosteroids.</li> <li>• People with poor co-ordination.</li> </ul>	

Key Features	
<b>Dexterity</b>	Low-moderate dexterity required
<b>Feedback: Dose Counter</b>	n/a
<b>Feedback: Taste/Feel</b>	May reduce taste of pMDI aerosol
<b>Feedback: Sound</b>	Click of valve audible when breathing in and out
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Consistent dosing when used correctly.
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Remove the cap from inhaler (and from the spacer if using an Aerochamber)
<b>2. Priming</b>	Shake the inhaler, then <b>place the mouthpiece of the inhaler into flat end of the spacer</b> and hold it upright, with your thumb on the base
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	<b>Place mouthpiece of spacer between teeth and lips and actuate one dose into spacer</b>
<b>5. Inhalation</b>	<b>Inhale slowly and deeply through the mouth</b> (if using Aero-chamber, a whistling sound indicates that the inspiratory rate is too fast)
<b>6. Breath holding</b>	Remove spacer from your mouth and hold breath for 10 seconds (or as long as comfortable), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Remove inhaler from the spacer and replace cap(s).


What you need to know / Critical Steps	
<b>Priming</b>	Delay between actuation of MDI into Spacer and inhalation results in reduced effective dose
<b>Multi-Dosing</b>	Multiple actuations of MDI into Spacer results in large particle sizes and reduced effective dose.
<b>Dose Wasting</b>	Important only to use single actuations into the spacer.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Slow and steady. Multiple breath
<b>Other information</b>	Useful for people with poor co-ordination. Useful when high doses of inhaled corticosteroid are prescribed
<b>Cleaning</b>	Check manufacturer's instructions. AeroChamber Plus: remove flat end of the spacer, and soak in lukewarm soapy water for a few minutes. Air dry in a vertical position, then reassemble when dry.

PODHALER	
	<b>Type of inhaler:</b> Single dose DPI - each dose contained in separate capsules
	<b>In use expiry:</b> Discard Podhaler device and case after 7 days of use
	<b>Available as</b>
<b>Antibiotic</b>	Tobi (tobramycin)

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	Remaining capsules may be counted.
<b>Feedback: Taste/Feel</b>	n/a
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	Capsules can be inspected to see if they are empty after inhalation.
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (40-85 L/min)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Wash and fully dry your hands. <b>Unscrew the top of the case. Then unscrew the mouthpiece.</b>
<b>2. Priming</b>	<b>Remove capsule from blister strip and insert capsule into inhaler chamber.</b> Screw the mouthpiece on until it stops. Do not overtighten. Hold the inhaler with the mouthpiece pointing down, then pierce the capsule by pressing the blue button <b>once and release.</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips, to make a tight seal</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 5 seconds, then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds, and breath normally, before <b>repeating steps 3-6 to ensure you inhaled the full dose.</b> Unscrew mouthpiece to empty the used capsule and take <b>the other 3 capsules in the same way.</b> When the full dose (4 capsules) have been inhaled, replace the mouthpiece and place the inhaler back in the storage case.

What you need to know / Critical Steps	
<b>Priming</b>	Capsules must only be removed immediately before use. Never place a capsule directly into the inhaler mouthpiece.
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Occasionally small capsule fragments may be inhaled (usually if the capsule is pierced more than once). This is not harmful.
<b>Cleaning</b>	Wipe the mouthpiece with a clean, dry cloth


PULVINAL	
	<b>Type of inhaler:</b> Reservoir Multidose DPI
	<b>In use expiry:</b> No restrictions
	<b>Available as</b>
<b>SABA</b>	Pulvinal salbutamol
<b>SAMA</b>	
<b>LABA</b>	
<b>LAMA</b>	
<b>LAMA/LABA</b>	
<b>ICS</b>	Pulvinal beclometasone
<b>ICS/LABA</b>	

Note: This inhaler device was discontinued in the UK in November 2015

Instructions for Use	
<b>1. Preparation</b>	<b>Unscrew cap.</b> Keeping the inhaler upright, gently tap it against a hard surface to level the powder inside the inhaler.
<b>2. Priming</b>	<b>Keeping the inhaler upright, press the button on the side of the inhaler and rotate the base as far as it will go in an anti-clockwise direction until a red mark shows in the hole below the button. Then whilst keeping the inhaler upright, twist the base back in a clockwise direction as far as it will go until it clicks and a green mark shows in the hole below the button.</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	Place mouthpiece between teeth and lips, whilst holding the inhaler upright.
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for a few seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use.

Key Features	
<b>Dexterity</b>	Moderate dexterity required
<b>Feedback: Dose Counter</b>	No dose counter, but medication powder is visible (replace when red ring at the bottom of the window is visible)
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (slow to fast)
<b>Device Resistance</b>	High airflow resistance
<b>Inspiratory Flow</b>	No validated test available

What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Risk of multi-dosing from multiple actuations
<b>Dose Wasting</b>	Exhaling into device will lose the dose
<b>Vents</b>	n/a
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Cap must be replaced after use to prevent exposure to moisture
<b>Cleaning</b>	Wipe the mouthpiece regularly with a soft cloth or fibre-free tissue

RESPIMAT		
	<b>Type of inhaler:</b> Soft Mist MDI	
	<b>In use expiry:</b> 3 months after canister has been loaded	
	<b>Available as</b>	
	SABA	Striverdi (olodaterol)
	SAMA	
	LAMA	Spiriva (tiotropium)
	LAMA/LABA	Spolto (tiotropium/olodaterol)
ICS		
ICS/LABA		

### Instructions for Use


#### Inserting the Cartridge

1. With cap closed, press safety catch and pull off clear base.
2. Push the cartridge into the inhaler completely & replace clear base.
3. Whilst holding the inhaler upright with the cap closed, turn the base in the direction of the red arrows until it clicks, and then open cap.
4. Point inhaler towards the ground, and press the dose release button.
5. Repeat steps 3-4 until a cloud is visible, then three further times (to fully prime the inhaler)

<b>1. Preparation</b>	Ensure cartridge has been inserted into the inhaler and the device primed for use correctly.
<b>2. Priming</b>	<b>Whilst holding the inhaler upright with the cap closed, turn the base in the direction of the red arrows until it clicks, and then open cap.</b>
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale slowly and deeply, pressing the button to release the medicine at the start of inhalation</b> , and continue to inhale deeply
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Discard capsule from central chamber and replace cap.

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	Has a dose indicator
<b>Feedback: Taste/Feel</b>	Taste of aerosol
<b>Feedback: Sound</b>	Audible 'click' on pressing the button to release a dose
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	Device locks when empty
<b>Dose Consistency</b>	Lung deposition may be reduced with fast inspiratory flow rates
<b>Device Resistance</b>	Low airflow resistance
<b>Inspiratory Flow</b>	No validated test available


What you need to know / Critical Steps	
<b>Priming</b>	Must be loaded and primed before first use. If not used for 7 days, spray one puff towards the ground; If not used for 21 days, re-prime device using spraying 3 puffs towards the ground
<b>Multi-Dosing</b>	Single dose - cannot multi-dose.
<b>Dose Wasting</b>	Ensure cap is closed when priming the device, to prevent inadvertent actuation of a dose
<b>Vents</b>	Do not block air vents
<b>Inspiratory Method</b>	Slow and steady
<b>Other information</b>	People with physical limitations may struggle to load the device
<b>Cleaning</b>	Clean the mouthpiece and metal part of the inhaler once a week with a damp cloth or tissue

SPIROMAX		
 <p>Photo. ©Teva</p>	<b>Type of inhaler:</b> Reservoir Multidose DPI	
	<b>In use expiry:</b> 6 months after removed from foil pouch	
	<b>Available as</b>	
	SABA	
	SAMA	
	LABA	
	LAMA	
LAMA/LABA		
ICS		
ICS/LABA	DuoResp (budesonide/formoterol)	

Key Features	
<b>Dexterity</b>	Low dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter. Counts down in steps of 2 doses.
<b>Feedback: Taste/Feel</b>	Taste of powder (contains lactose), and gritty texture of powder
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Increased dose delivery at higher inspiratory flow rates (90L/min vs. 40L/min), but likely insignificant clinical impact.
<b>Device Resistance</b>	Medium airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use.
<b>2. Priming</b>	Hold inhaler in the horizontal or upright position, and shake the device. <b>Open the cap completely.</b>
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips
<b>5. Inhalation</b>	Inhale strongly and deeply.
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Close cover after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the horizontal or vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose. Multiple actuations of device winds dose counter down, but does not load dose (appears empty before it actually is).
<b>Dose Wasting</b>	Exhaling into device will lose the dose
<b>Vents</b>	Do not block air vents
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Dose emission from Spiromax is less dependent on airflow than the Turbohaler
<b>Cleaning</b>	Wipe the mouthpiece with a dry tissue

TURBOHALER		
 <p>Photo. ©AstraZeneca</p>	<b>Type of inhaler:</b> Reservoir Multidose DPI	
	<b>In use expiry:</b> No restrictions	
	<b>Available as</b>	
	SABA	Bricanyl (terbutaline)
	SAMA	
	LABA	Oxis (formoterol)
	LAMA	
	LAMA/LABA	
ICS	Pulmicort (budesonide)	
ICS/LABA	Symbicort (budesonide/formoterol)	

Key Features	
<b>Dexterity</b>	Moderate dexterity required
<b>Feedback: Dose Counter</b>	Has a dose counter (Symbicort), which counts in steps of 20 doses; or dose indicator (Bricanyl, Oxis, Pulmicort), indicates empty when red mark reaches bottom of indicator window.
<b>Feedback: Taste/Feel</b>	Generally no taste
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Inconsistent dosing across range of inspiratory flow rates (lower effective dose with low inhalation speed)
<b>Device Resistance</b>	Medium airflow resistance (Symbicort); Medium-high airflow resistance (Pulmicort, Oxis, Bricanyl)
<b>Inspiratory Flow</b>	Test with In Check DIAL inspiratory flow meter, or with Turbohaler whistle.

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use, and then unscrew cap.
<b>2. Priming</b>	Shake the inhaler, then whilst <b>keeping it upright, rotate grip as far as it will go in one direction. Then turn it as far as it will go in the other direction</b> (it does not matter which way you turn it first). You should hear a click sound.
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use.


What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position
<b>Multi-Dosing</b>	Single dose - cannot multi-dose. Multiple actuations of device winds dose counter down, but does not load dose (appears empty before it actually is)
<b>Dose Wasting</b>	Exhaling into device will lose the dose
<b>Vents</b>	Do not block air vents
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Turn Aid available for people with physical limitations. Desiccant can be heard when shaking the device - it is not the amount of drug remaining.
<b>Cleaning</b>	Wipe the mouthpiece with a dry tissue once a week

TURBOSPIN	
 Photo. ©Actavis	<b>Type of inhaler:</b> Single dose DPI - each dose contained in separate capsules
	<b>In use expiry:</b> Discard Turbospin inhaler after completion of each treatment pack
<b>Available as</b>	
<b>Antibiotic</b>	Colobreathe (colistimethate sodium)

Instructions for Use	
<b>1. Preparation</b>	Remove the cap. Then unscrew the mouthpiece.
<b>2. Priming</b>	Remove capsule from blister strip and insert capsule into inhaler chamber with the widest end first. Screw the mouthpiece on until it stops. Do not overtighten. Hold the inhaler with the mouthpiece pointing upwards then pierce the capsule by gently pushing the piston upwards <b>once and release.</b>
<b>3. Exhaling</b>	Exhale fully and away from mouthpiece
<b>4. Mouth</b>	Place mouthpiece between teeth and lips, to make a tight seal
<b>5. Inhalation</b>	Inhale slowly and deeply, but at a rate sufficient to hear the capsule vibrate. (If the capsule does not vibrate, the chamber of the inhaler device should be tapped gently and the inhalation repeated.)
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for about 10 seconds, then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds, and breath normally, before repeating steps 3-6 to ensure you inhaled the full dose. Unscrew mouthpiece to empty the used capsule then replace the mouthpiece and replace the cap.

Key Features	
<b>Dexterity</b>	High dexterity required
<b>Feedback: Dose Counter</b>	Remaining capsules may be counted.
<b>Feedback: Taste/Feel</b>	Taste of powder (may be unpleasant)
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	Capsules can be inspected to see if they are empty after inhalation.
<b>Device Lock-Out</b>	n/a
<b>Dose Consistency</b>	Relatively consistent dosing across range of inspiratory flow rates (30-90 L/min)
<b>Device Resistance</b>	Medium-low airflow resistance
<b>Inspiratory Flow</b>	No validated test available

What you need to know / Critical Steps	
<b>Priming</b>	Capsules must only be removed immediately before use. Do not pierce the capsule more than once
<b>Multi-Dosing</b>	Single dose - cannot multi-dose
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	Do not cover the air slits with your fingers or mouth during inhalation.
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Occasionally small capsule fragments may be inhaled (usually if the capsule is pierced more than once). This is not harmful.
<b>Cleaning</b>	Clean after each dose. Press the piston down firmly a few times, whilst keeping the chamber upside down. Clean the chamber using a tissue or cotton bud, then screw the mouthpiece back on and replace the cap.

TWISTHALER	
	<b>Type of inhaler:</b> Reservoir Multidose DPI
	<b>In use expiry:</b> 3 months after removed from foil pouch
	<b>Available as</b>
SABA	
SAMA	
LABA	
LAMA	
LAMA/LABA	
ICS	Asmanex (mometasone)
ICS/LABA	

Key Features	
<b>Dexterity</b>	Low-moderate dexterity required
<b>Feedback: Dose Counter</b>	Dose counter. Counts individual doses. Shaded indicator shows for last 20 doses.
<b>Feedback: Taste/Feel</b>	Generally no taste
<b>Feedback: Sound</b>	n/a
<b>Feedback: Visual</b>	n/a
<b>Device Lock-Out</b>	Cap locks in place when empty
<b>Dose Consistency</b>	Consistent dosing across range of inspiratory flow rates (28-70 L/min)
<b>Device Resistance</b>	Medium-high airflow resistance
<b>Inspiratory Flow</b>	No validated test available

Instructions for Use	
<b>1. Preparation</b>	Check dose counter before use, and ensure that the counter and the pointer on the cap are lined up.
<b>2. Priming</b>	Shake the inhaler, then whilst <b>keeping it upright, grip the base and twist the cap anticlockwise to remove it. Ensure that the counter and the pointer on the inhaler device are now lined up.</b>
<b>3. Exhaling</b>	<b>Exhale fully and away from mouthpiece</b>
<b>4. Mouth</b>	<b>Place mouthpiece between teeth and lips</b>
<b>5. Inhalation</b>	<b>Inhale strongly and deeply.</b>
<b>6. Breath holding</b>	Remove inhaler from your mouth and hold your breath for up to 10 seconds (or as long as possible), then breathe out slowly.
<b>7. Closing &amp; Repeating</b>	Wait a few seconds before repeating steps 3-6 to ensure you inhaled the full dose. Replace cap after use, and ensure that the counter and the pointer on the cap are lined up again.

What you need to know / Critical Steps	
<b>Priming</b>	Must be primed in the vertical position. Removal of cap primes the device
<b>Multi-Dosing</b>	Multiple actuations of device winds dose counter down, but does not load dose (appears empty before it actually is)
<b>Dose Wasting</b>	Exhaling into device will lose the dose. Tilting, inverting or shaking the device may lose the dose.
<b>Vents</b>	Do not block the air vents
<b>Inspiratory Method</b>	Quick and deep
<b>Other information</b>	Removing the cap causes the dose counter to count down
<b>Cleaning</b>	Wipe the outside of the mouthpiece with a dry cloth or tissue.

## Monograph Bibliography

1. Anonh A, Leinonen M et Ranki-Pesonen M. Patient satisfaction with Easyhaler® compared with other inhalation systems in the treatment of asthma: a meta-analysis. *Current Therapeutic Research* 2000;61:61-73.
2. Al-Showair RAM, Pearson SB et Chrystyn H. The potential of a 2Tone Trainer to help patients use their metered-dose inhalers. *Chest* 2007;131:1776-82.
3. Al-Showair RAM, Tarsin WJ, Assi KH et al. Can all patients with COPD use the correct inhalation flow with all inhalers and does training help? *Respir Med* 2007;101:2395-401.
4. Amnari WG et Chrystyn H. Optimizing the Inhalation Flow and Technique Through Metered Dose Inhalers of Asthmatic Adults and Children Attending a Community Pharmacy. *Journal of Asthma* 2013;50:505-13.
5. Anderson P. Patient preference for and satisfaction with inhaler devices. *European Respiratory Review*. 2005;14:109-16.
6. Asthma UK. Using your inhalers. Last reviewed June 2015. Available at: [www.asthma.org.uk/advice-using-inhalers](http://www.asthma.org.uk/advice-using-inhalers)
7. Azouz W et Chrystyn H. Clarifying the dilemmas about inhalation techniques for dry powder inhalers: integrating science with clinical practice. *Primary Care Respiratory Journal* 2012;21:208-13.
8. Azouz W, Chetcuti P, Hosker et al. Inhalation characteristics of asthma patients, COPD patients and healthy volunteers with the Spiromax® and Turbuhaler® devices: a randomized, cross-over study. *BMC Pulmonary Medicine* 2015;15:4.
9. Barry PW et al. Inhalational drug delivery from seven different spacer devices. *Thorax* 1996;51:835-40.
10. Barry PW et O'Callaghan C. Multiple actuations of salbutamol MDI into a spacer device reduce the amount of drug recovered in the respirable range. *Eur Respir J* 1994;7:1707-9.
11. Bashedi IA, Armour CL, Bosnic-Anticevich SZ et al. Evaluation of a novel educational strategy, including inhaler-based reminder labels, to improve asthma inhaler technique. *Patient Education and Counseling* 2008;72:16-33.
12. Bashedi IA, Bosnic-Anticevich SZ, Armour CL et al. Checklists for Powder Inhaler Technique: A Review and Recommendations. *Respiratory Care* 2014;59:140-54.
13. Bashedi IA, Bosnic-Anticevich SZ et al. User error with Diskus and Turbuhaler by asthma patients and pharmacists in Jordan and Australia. *Respiratory Care* 2013;58:1916-23.
14. Bashedi IA, Reddel HK, Armour CL et al. Counseling about turbuhaler technique: needs assessment and effective strategies for community pharmacists. *Respir Care* 2005;50:617-23.
15. Bashedi IA, Reddel HK, Armour CL et al. Improved asthma outcomes with a simple inhaler technique intervention by community pharmacists. *Journal of Allergy & Clinical Immunology* 2007;119:1537-8.
16. Baverstock M, Woodhall N et Maarman V. P94 Do healthcare professionals have sufficient knowledge of inhaler techniques in order to educate their patients effectively in their use. *Journal of Clinical Pharmacy and Therapeutics* 2017;42:17-18.
17. Bell J. Why are inhaler techniques in asthma and COPD? *The British Journal of Primary Care Nursing*. 2008;2:37-39.
18. Borgström L, Bondesson E, Morén F et al. Lung deposition of budesonide inhaled via Turbuhaler®: a comparison with terbutaline sulphate in normal subjects. *Eur Respir J* 1994;7:69-73.
19. Brand P, Hedere B, Austen G et al. Higher lung deposition from RespiMat/Soft Mist Inhaler (SMI) than HFA-MDI in COPD patients with poor technique. *Journal of Aerosol Medicine* 2007;20:165.
20. Brand P, Hedere B, Austen G et al. Higher lung deposition with RespiMat® Soft Mist™ Inhaler than HFA-MDI in COPD patients with poor technique. *Int J Chron Obstruct Pulmon Dis* 2008;3(4):763-70.
21. Brennan VK, Osman LM, Graham H et al. True device compliance: the need to consider both competence and contrivance. *Respir Med* 2005;99:97-102.
22. Brocklebank D, Ram F, Wright J et al. Comparison of the effectiveness of inhaler devices in asthma and chronic obstructive airways disease: a systematic review of the literature. *Health Technology Assessment* 2001;5:1-155.
23. Broeders ME, Vincken W, Corbetta L. The ADMIT series—issues in inhalation therapy. 7. Ways to improve pharmacological management of COPD: the importance of inhaler choice and inhalation technique. *Primary Care Respiratory Journal* 2011;20:338-343.
24. Brown PH, Blundell G, Greening AP et al. Do large volume spacer devices reduce the systemic effects of high dose inhaled corticosteroids? *Thorax* 1990;45:736-9.
25. Caillaud C, Le Merre C, Martinat Y et al. Dose-dependent response to tiotropium administered via the RespiMat® Soft Mist™ Inhaler (SMI) to patient with COPD. *European Respiratory Journal* 2006;28(suppl 50):318S-9S.
26. Canonica GW, Arp J, Keegstra JR et al. Spiromax, a new dry powder inhaler: dose consistency under simulated real-world conditions. *Journal of Aerosol Medicine and Pulmonary Drug Delivery* 2015;28:309-19.
27. Capstick TGD et Clifton IJ. Inhaler technique and training in people with chronic obstructive pulmonary disease and asthma. *Expert Rev Respir Med* 2012;6:91-103.
28. Chapman KR, Fogarty CM, Peckitt C et al. Delivery characteristics and patients' handling of two single-dose dry powder inhalers used in COPD. *International Journal of COPD* 2011;6:365-63.
29. Chodosh S, Flanders JS, Kesten S et al. Effective Delivery of Particles with the HandiHaler® Dry Powder Inhalation System over a Range of Chronic Obstructive Pulmonary Disease Severity. *J Aerosol Med* 2001;14:309-15.
30. Chrystyn H et Niederlaender C. The Genuair® inhaler: a novel, multidose dry powder inhaler. *Int J Clin Pract* 2012;66:309-17.
31. Chrystyn H et Price D. Not all asthma inhalers are the same: factors to consider when prescribing an inhaler. *Prim Care Respir J* 2009;18:243-9.
32. Chrystyn H. Is inhalation rate important for a dry powder inhaler? using the In-Check Dial to identify these rates. *Respir Med* 2003;97:181-7.
33. Cochran JJ, Boudry MV, Downs KE et al. Inhaled Corticosteroids for Asthma Therapy: Patient Compliance, Devices, and Inhalation Technique. *Chest* 2000;117:542-50.
34. Conole D et Keating GM. colistimethate sodium drug powder for inhalation: a review of its use in the treatment of chronic Pseudomonas aeruginosa infection in patients with cystic fibrosis. *Drugs* 2014;74:377-87.
35. Corradi M, Chrystyn H, Cosio BG et al. NEXThaler, an innovative dry powder inhaler delivering an extrafine fixed combination of beclomethasone and formoterol to treat large and small airways in asthma. *Expert Opin Drug Deliv*. 2014;11:1497-1506.
36. Crompton GK, Barnes PJ, Broeders M et al. The need to improve inhalation technique in Europe: a report from the Aerosol Drug Management Improvement Team. *Respir Med* 2009;103:79-94.
37. Dahl R, Backer V, Allgaard B et al. Assessment of patient preference of the handihaler compared with the metered dose inhaler four weeks after instruction. *Respiratory Medicine* 2003;97:1126-33.
38. Dal Negro R, Micheletto C, Tognella S et al. Evidence of Adequacy of the Performance of the Pulvinal™ by Measuring Through-Device Peak Inspiratory Flow Rate in Severe Airways Obstruction in Adults and Children. *J Aerosol Med* 2001;14:343-9.
39. Demoly P, Hagedoorn P, de Boer AH et al. The clinical relevance of dry powder inhaler performance for drug delivery. *Respir Med* 2014;108:1195-1203.
40. Dolovich MB, Ahrens RC, Hess DR et al. Device selection and outcomes of aerosol therapy: Evidence-based guidelines: American College of Chest Physicians/American College of Asthma, Allergy, and Immunology. *Chest* 2005;127:335-71.
41. Doyle S, Lloyd A, Williams A et al. What happens to patients who have their asthma device switched without their consent? *Primary Care Respiratory Journal* 2010;19:131-9.
42. Erickson SR, Landino HM, Zarowitz BJ et Kirkl DM. Pharmacists' understanding of patient education on metered-dose inhaler technique. *The Annals of Pharmacotherapy* 2000;34, No. 11, pp. 1249-1256
43. European Medicines Agency. Committee for Medicinal Products for Human Use (CHMP). CHMP assessment report: DuoResp Spiromax. EMA/CHMP/17092/2014. Available at: [http://www.ema.europa.eu/docs/en\\_GB/document\\_library/EPAR\\_-\\_Public\\_assessment\\_report/human/000234/WC500167183.pdf](http://www.ema.europa.eu/docs/en_GB/document_library/EPAR_-_Public_assessment_report/human/000234/WC500167183.pdf). Accessed 25th February 2014.
44. Everard M. Guidelines for devices and choices. *J Aerosol Med* 2001;14 Suppl 1:S59-S64.
45. Everard ML, Devadason SG, Le Souëf PN. Flow early in the inspiratory manoeuvre affects the aerosol particle size distribution from a Turbuhaler. *Respir Med* 1997;91:624-8.
46. Fenton C, Keating GM et Plosker GL. Novolizer®: A Multidose Dry Powder Inhaler. *Drugs* 2003;63:2437-45.
47. Geller DE, Weers J et Heuveling S. development of an inhaled dry-powder formulation of tobramycin using PulmoSphere™ technology. *Journal of Aerosol Medicine and Pulmonary Drug Delivery* 2011;24:175-82.
48. Giraud V et al. Prime D. In vivo product performance of fluticasone furoate/vilanterol delivered from a novel dry powder inhaler, using the electronic lung to replicate asthma and COPD patient inhalation profiles. *Am J Respir Crit Care Med*. 2012;185:A2940.
49. Giraud V et Roche N. Misuse of corticosteroid metered-dose inhaler is associated with decreased asthma stability. *European Respiratory Journal* 2002;19:246-51.
50. GlaxoSmithKline. Multiple indexing of Relvar Ellipta without inhalation. Personal Communication. 11th May 2015.
51. Global Initiative for Asthma. Instructions for Inhaler and Spacer Use. Available at: <http://www.ginasthma.org/documents/11>
52. Gustafsson P, Taylor A, Zanen P et al. Can patients use all dry powder inhalers equally well? *Int J Clin Pract* 2005;59(Suppl. 149):13-8.
53. Hamilton J, Prime D et Haggret R. Ex vivo product performance of fluticasone furoate/vilanterol delivered from a novel dry powder inhaler, using the electronic lung to replicate asthma and COPD patient inhalation profiles. *Am J Respir Crit Care Med*. 2012;185:A2940.
54. Hanania NA, Wittman R, Kesten S et al. Medical personnel's knowledge of and ability to use inhaling devices: metered-dose inhalers, spacing chambers and breath-actuated dry powder inhalers. *Chest* 1994;105:111-7.
55. Haughney J, Price D, Barnes NC, et al. Choosing inhaler devices for people with asthma: current knowledge and outstanding research needs. *Respiratory Medicine*. 2010;104:1237-45.
56. Hindle M et Chrystyn H. Relative bioavailability of salbutamol to the lung following inhalation using metered dose inhalation methods and spacer devices. *Drugs* 1994;48:53-53.
57. Hirst PH, Bacon RE, Pitcairn GR et al. A comparison of the lung deposition of budesonide from Easyhaler®, Turbuhaler® and pMDI plus spacer in asthmatic patients. *Respiratory Medicine* 2001;95:720-7.
58. Janssens W, VandenBrande P, Hardeman E et al. Inspiratory flow rates at different levels of resistance in elderly COPD patients. *Eur Respir J* 2008;31:78-83.
59. Kanabuchi K, Kondo T, Tanigaki T et al. Minimal inspiratory flow from dry powder inhalers according to a biphasic model of pressure vs. flow relationship. *Tokai J Exp Clin Med* 2011;36:1-4.
60. Killefether S, Pontiz H, Beck E et al. Improved delivery of ipratropium bromide/fenoterol from RespiMat® Soft Mist™ inhaler in patients with COPD. *Respiratory Medicine* 2004;98:387-97.
61. Koskela T, Malmstrom K, Sairanen et al. Efficacy of salbutamol via Easyhaler® unaffected by low inspiratory flow. *Respiratory Medicine* 2000;94:1229-33.
62. Kuna P, Thyroff-Friesinger U, Gath I et al. Randomized equivalence trial: A novel multidose dry powder inhaler and originator device in adult and adolescent asthma. *Allergy Asthma Proc* 2015;36:352-64.
63. Lavorini F, Levy ML, Corrigan C et al. The ADMIT series – issues in inhalation therapy. 6) Training tools for inhalation devices. *Primary Care Respiratory Journal* 2010;19:335-41.
64. O'Callaghan C, Magan A, Dubus JC et al. Effect of incorrect use of dry powder inhalers on management of patients with asthma and COPD. *Respir Med* 2008;102:593-604.
65. Lenney J, Innes J, Crompton G. Inappropriate inhaler use: assessment of use and patient preference of seven inhalation devices. *Respiratory Medicine*. 2000;94:496-500.
66. Levy ML, Hardwell A, McKnight E et al. Asthma patients' inability to use a pressurised metered-dose inhaler (pMDI) correctly correlates with poor asthma control as defined by the Global Initiative for Asthma (GINA) strategy: a retrospective analysis. *Primary Care Respiratory Journal* 2013;22:406-11.
67. Magnusson H, Watz H, Zimmermann N et al. Peak inspiratory flow through the Genuair® inhaler in patients with moderate or severe COPD. *Respir Med* 2009;103:1832-7.
68. Malmstrom K, Sorva R, Silvestri M et al. Application and efficacy of the multi-dose powder inhaler, Easyhaler®, in children with asthma. *Pediatric Allergy and Immunology* 1999;10:66-70.
69. McKee DC. Tobramycin Inhalation Powder: A review of its use in the treatment of chronic Pseudomonas aeruginosa infection in patients with cystic fibrosis. *Drugs* 2013;73:1815-37.
70. Melani AS, Bonavia M, Clienti V et al. Inhaler mishandling remains common in real life and is associated with reduced disease control. *Respir Med* 2011;105:930-8.
71. Melani AS, Bonavia M, Melletta D, Barbato N et al. Inhalation technique and variables associated with misuse of conventional metered-dose inhalers and newer dry powder inhalers in experienced adults. *Annals of Allergy, Asthma, & Immunology* 2004;93:439-46.
72. Molimard M, Raheerich C, Lignot S et al. Assessment of handling of inhaler devices in real life: an observational study in 3811 patients in primary care. *J Aerosol Med* 2003;16:249-254.
73. Newhouse MT, Nantel NP, Chambers CB et al. Clickhaler (a Novel Dry Powder Inhaler) Provides Similar Bronchodilation to Pressurized Metered-Dose Inhaler, Even at Low Flow Rates. *Chest* 1999;115:952-5.
74. Newman S, Stred K, Hooper G et al. Comparison of gamma scintigraphy and a pharmacokinetic technique for assessing pulmonary deposition of terbutaline sulphate delivered by pressurized metered dose inhaler. *Pharmaceutical Research* 1995;12:231-6.
75. Newman SP, Pitcairn GR, Hirst PH et al. Scintigraphic comparison of budesonide deposition from two dry powder inhalers. *Eur Respir J* 2000;16:178-183.
76. Newman SP, Pitcairn GR, Adkin DA et al. Comparison of beclomethasone dipropionate delivery by Easyhaler® dry powder inhaler and pMDI plus large volume spacer. *Journal of Aerosol Medicine* 2001;14: 217-25.
77. Newman SP, Weisz AWB, Talaei N et al. Improvement of drug delivery with a breath actuated pressurised aerosol for patients with poor inhaler technique. *Thorax* 1991;46:712-6.
78. Newman SP. Aerosol deposition considerations in inhalation therapy. *Chest* 1985;88:1525-160S.
79. Pitcairn G, Reader S, Pavia D et Newman S. Deposition of corticosteroids aerosol in the human lung by RespiMat® Soft Mist™ inhaler compared to deposition by metered dose inhaler or by Turbuhaler® dry powder inhaler. *Journal of Aerosol Medicine* 2005;18:264-272.
80. Nsour WM, Aldred A, Corrado OJ et al. Measurement of peak inhalation rates with an In-Check Meter® to identify an elderly patient's ability to use a Turbuhaler®. *Respir Med* 2001;95:965-8.
81. O'Callaghan C, Cant M et Robertson C. Delivery of beclomethasone dipropionate from a spacer device: what dose is available for inhalation? *Thorax* 1994;49:961-4.
82. Palander A, Mattila T, Karhu M et al. In vitro comparison of three salbutamol-containing multidose dry powder inhalers. *Clinical Drug Investigation* 2000;20:25-33.
83. Pasquali I, Brambilla G et Copelli D. Effect of flow rate on dose delivery of three dry powder inhalers: NEXThaler®, Turbuhaler® and Diskus®. *Respiratory Drug Delivery Europe* 2013;2:267-72.
84. Pavkov R, Mueller S, Fiebach K et al. Characteristics of a capsule based dry powder inhaler for the delivery of inhaled corticosteroid. *Current Medical Research & Opinion* 2010;26:257-63.
85. PHT® Turbospin™ dry powder inhaler. Brochure. Available at: <http://www.phtpharma.com/medical-devices/products/turbospin-dpi> Accessed 20th October 2015.
86. Pitcairn G, Lunghetti G, Ventura P et al. A comparison of the lung deposition of salbutamol inhaled from a new dry powder inhaler, at two inhaled flow rates. *International journal of pharmaceuticals* 1994;102:11-18.
87. Pitcairn G, Reader S, Pavia D et Newman S. Deposition of corticosteroids aerosol in the human lung by RespiMat® Soft Mist™ inhaler compared to deposition by metered dose inhaler or by Turbuhaler® dry powder inhaler. *Journal of Aerosol Medicine* 2005;18:264-272.
88. Price D, Bosnic-Anticevich S, Briggs A et al. Inhaler competence in asthma: Common errors, barriers to use and recommended solutions. *Respiratory Medicine* 2013;107:37-46.
89. Price D, Thomas WJ, Mitchell G et al. Improvement of asthma control with a breath-actuated pressurised metered dose inhaler (BAI): a prescribing claims study of 556 patients using a traditional pressurised metered dose inhaler (MDI) or a breath-actuated metered dose inhaler. *Respir Med* 2003;97:12-9.
90. Prihe D, De Backer W, Hamilton M et al. Comparison of inhalation profiles through a novel dry powder inhaler (ndpi) and lung function measurements for healthy subjects, asthma and chronic obstructive pulmonary disease (COPD) patients. *Am J Respir Crit Care Med*. 2012;185:A2941.
91. Riley J, Tabberer M, Richard N et al. Use of a new dry powder inhaler to deliver inhaled corticosteroid in the treatment of COPD. *European Respiratory Society Congress 2013*. Poster No. P4145. Available at: [http://hrs.ersnet.org/abstract\\_print\\_13/files/395.pdf](http://hrs.ersnet.org/abstract_print_13/files/395.pdf)
92. Rönmark C, Erik R, Lindqvist A et al. Correct use of three powder inhalers: comparison between diskus, turbuhaler, and easyhaler. *J Asthma* 2005;42:173-8.
93. Royal College of Physicians. Why asthma still kills: the national review of asthma deaths (NRAD) confidential enquiry report. London: RCP. 2014.
94. Schulte M, Osseiran K, Betz R et al. Handling of and preferences for available dry powder inhaler systems by patients with asthma and COPD. *Journal of Aerosol Medicine and Pulmonary Drug Delivery* 2008;21:321-8.
95. Schultz A, Schulte J, Venter A et al. Aerosol inhalation from spacers and valved holding chambers requires few tidal breaths for children. *Pediatrics* 2010;126:e1493-e1498.
96. Schürmann W, Eichmüdner M, Moroni P et al. RespiMat® Soft Mist™ Inhaler versus hydrofluoroalkane metered dose inhaler: patient preference and satisfaction. *Treatments in Respiratory Medicine* 2005;4:53-61.
97. Schuster A, Halliburton C, Döring G et al. Safety, efficacy and convenience of colistimethate sodium dry powder for inhalation (Colobreathe DPI) in patients with cystic fibrosis: a randomised study. *Thorax* 2013;68:344-50.
98. Schweisfurth H, Mallinen A, Koskela T et al. Comparison of two budesonide powder inhalers, Easyhaler® and Turbuhaler®, in steroid-naïve asthmatic patients. *Respiratory Medicine* 2002;96:599-606.
99. Simple Steps Education: Inhaler Technique – 7 Steps to Success Reminder Cards. [www.simplestepseducation.co.uk](http://www.simplestepseducation.co.uk)
100. Singh D, Fiebach K, Dederichs J et al. Dose delivery characterization of inhaled corticosteroid delivered from diskus following inhalation by COPD patients. *Am J Respir Crit Care Med* 2010;181:A4419.
101. Summary of Product Characteristics, Electronic Medicines Compendium. Datapharm Communications Ltd. <http://emc.medicines.org.uk/102>
102. Sveasater H, Dale P, Garrill K et al. Qualitative assessment of attributes and ease of use of the ELLIPTA™ dry powder inhaler for delivery of maintenance therapy for asthma and COPD. *BMC Pulmonary Medicine* 2013;13:72.
103. Taylor A et Gustafsson P. Do all dry powder inhalers show the same pharmaceutical performance? *Int J Clin Pract* 2005;59(Suppl. 149):7-12.
104. Thomas WJ, Price O, Chrystyn H et al. Inhaled corticosteroids for asthma: impact of practice level device switching on asthma control. *BMC Pulmonary Medicine* 2009;9:9.
105. Thompson J, Irvine T, Grathwohl K et Roth B. Misuse of metered-dose inhalers in hospitalized patients. *Chest* 1994;105:715-717.
106. Tukiainen H, Ryttila P, Hamalainen KM et al. Safety, tolerability and acceptability of two dry powder inhalers in the administration of budesonide in steroid-treated asthmatic patients. *Respiratory Medicine* 2002;96:221-9.
107. van der Palen J, Ginko T, Kroker A et al. Preference, satisfaction and errors with two dry powder inhalers in patients with COPD. *Expert Opin. Drug Deliv* 2013;10:1623-31.
108. van der Palen J, Klein JJ, Schildkamp AM. Comparison of a new multi-dose powder inhaler (Diskus/Accuhaler) and the Turbuhaler regarding preference and ease of use. *J Asthma* 1998;35:147-52.
109. van der Palen J. Genuair® in chronic obstructive pulmonary disease: a novel, user-friendly, multidose, dry powder inhaler. *Ther Deliv* 2014;5:795-800.
110. Van Noord JA, Cornelissen PJG, Aumann, JL et al. Efficacy in COPD patients of tiotropium administered with the RespiMat® Soft Mist™ inhaler compared to Handihaler®. *European Respiratory Journal* 2006;28(suppl 50):431S.
111. VanDevanter DR et Geller DE. Tobramycin administered by the TOBI® Podhaler® for persons with cystic fibrosis: a review. *Medical Devices: Evidence and Research* 2011;4:179-188.
112. Vanto T, Hamalainen KM, Vaiteristo M et al. Comparison of two budesonide dry powder inhalers in the treatment of asthma in children. *Journal of Aerosol Medicine* 2004;17:15-24.
113. Vidgren M, Silvestri M, Vidgren P et al. Easyhaler® multiple dose powder inhaler - practical ad effective alternative to the pressurised MDI. *Aerosol Science and Technology* 1995;22:335-45.
114. Vincken W, Bantje T, Middle M et al. Long-term efficacy and safety of ipratropium bromide plus fenoterol via RespiMat® Soft Mist™ inhaler (SMI) versus a pressurised metered-dose inhaler in asthma. *Clinical Drug Investigation* 2004;24:7-28.
115. Virchow JC, Crompton GK, Dal Negro et al. Importance of inhaler devices in the management of airway disease. *Respir Med* 2008;102:10-9.
116. Virchow JC, Weuthen T, Harmer OJ et al. Identifying the features of an easy-to-use and intuitive dry powder inhaler for asthma and chronic obstructive pulmonary disease therapy: results from a 28-day device handling study, and an airflow resistance study. *Expert Opinion on Drug Delivery* 2014;11:1849-57.
117. Virchow JC, Palani M, Linnane P et al. Comparing usability of NEXThaler® with other inhaled corticosteroid/long-acting beta2-agonist fixed combination dry powder inhalers in asthma patients. *Journal of aerosol medicine and pulmonary drug delivery* 2014;27:1-8.
118. Yang TT, Li S, Wyka B et al. Drug delivery performance of the mometasone furoate dry powder inhaler. *J Aerosol Med* 2001;14:487-94.
119. ADMIT. Inhalers. Available at: <http://www.admit-inhalers.org/en/>
120. MHRA. Medicines Information: SPC & PILs. Available at: <http://www.mhra.gov.uk/spc/pil/>